

## Preparation of complete medium for malaria culture v1.2

Procedure



***In vitro* Module**

**WorldWide Antimalarial Resistance Network (WWARN)**



**Suggested citation:** *In vitro* Module, WWARN. 2011. Preparation of complete medium for malaria culture. WWARN procedure.

**Procedure ID:** INV02

**This procedure was developed by:**

*In Vitro* Module, WWARN

#### Version History

Version number	Revision(s) & reason for amendment	Release Date
1.0	Creation of the document	16/08/2010
1.1	Changes to Sections 6.2 and 6.3	01/04/2011
1.2	Changes to template	31/08/2012

**For more information, contact:**

invitro@wwarn.org

WorldWide Antimalarial Resistance Network (WWARN)

[www.wwarn.org](http://www.wwarn.org)

## Contents

1. Purpose .....	4
2. Scope .....	4
3. Abbreviations .....	4
4. Duties and Responsibilities .....	4
5. Materials and Equipment .....	4
5.1 Materials .....	4
5.2 Equipment .....	5
6. Procedures .....	5
6.1 5% sodium bicarbonate solution .....	5
6.2 RPMI washing medium .....	5
6.3 Complete medium .....	6
6.4 Quality control (QC) .....	6
7. References .....	7

## 1. Purpose

This protocol describes the preparation of medium used to culture *P. falciparum* isolates or clones.

## 2. Scope

This procedure is designed for use by laboratories testing *in vitro* drug susceptibility using *P. falciparum*. This procedure is applicable in well-equipped cell culture laboratories. Training is required to perform the procedure successfully.

## 3. Abbreviations

RPMI	Roswell Park Memorial Institute medium 1640
RT	Room temperature
FW	Formula weight

## 4. Duties and Responsibilities

The preparation of complete medium for malaria cultures must be carried out by a competent technician.

## 5. Materials and Equipment

### 5.1 Materials

- Sterile graduated pipettes
- 100mL volumetric flask
- 1L volumetric flask
- Cups and spatula
- Disposable 0.22µM sterilisation filter unit and filter unit receiver
- Sterile vials
- RPMI 1640, powdered medium, stored at 4°C
- Sodium bicarbonate, NaHCO<sub>3</sub>, FW 84.01, stored at RT
- HEPES, C<sub>8</sub>H<sub>18</sub>N<sub>2</sub>O<sub>4</sub>S, FW 238.1, stored at RT
- Hypoxanthine, C<sub>5</sub>H<sub>4</sub>N<sub>4</sub>O, FW 136.11, stored at RT
- L-glutamine, C<sub>5</sub>H<sub>10</sub>N<sub>2</sub>O<sub>3</sub>, FW 146.1, stored at RT
- Glucose, C<sub>6</sub>H<sub>12</sub>O<sub>6</sub>, FW 180.16, stored at RT
- Orotic acid, C<sub>5</sub>H<sub>4</sub>N<sub>2</sub>O<sub>4</sub>, FW156.1 (optional), stored at RT
- Gentamicin (optional), C<sub>21</sub>H<sub>43</sub>N<sub>5</sub>O<sub>7</sub>, FW 477.596, stored at -20°C

- Sterile human serum stored at -20 °C or Albumax stored at + 4 °C, protected from extended exposure to light
- Hydrochloric acid 1N
- Sodium hydroxide 1N
- Water for cell culture or double-distilled water

## 5.2 Equipment

- Analytical balance
- Magnetic stir plate
- Laminar flow hood
- Vacuum trap
- Cryogenic equipment at + 4 °C
- Cryogenic equipment at – 20 °C
- Water-bath or heater block
- pH meter
- Osmometer

## 6. Procedures

### 6.1 5% sodium bicarbonate solution

- I. Weigh 5g sodium bicarbonate
- II. Add 100mL of cell culture water in volumetric flask
- III. Stir on magnetic stirring plate until dissolved
- IV. Store at room temperature
- V. Use within 24 hours of preparation

### 6.2 RPMI washing medium

Used to wash samples and red blood cells.

- I. Weigh 10.43g of powdered RPMI medium and add to 1L volumetric flask
- II. Add 5.95g HEPES (25mM final concentration)
- III. Add 42mL 5% sodium bicarbonate (as prepared in 6.1) (25mM final concentration)
- IV. Fill to 1L with cell culture water
- V. Stir the solution until dissolved
- VI. Using a calibrated meter, check medium pH. If necessary, adjust to  $\text{pH}7.0 \pm 0.3$  with either hydrochloric acid (1N) or sodium hydroxide (1N). The medium should be orange, not red/purple or yellow, in colour due to the phenol red dye
- VII. Using a calibrated meter, check that the medium osmolarity =  $292 \pm 5\%$  mOsm. If necessary, adjust with sterile water
- VIII. In a sterile environment, use a 0.22 $\mu\text{M}$  filter to sterilise the medium
- IX. Aseptically dispense medium into sterile bottles
- X. Store at – 20 °C for a maximum 6 months
- XI. After first use, store at + 4 °C for one week

### 6.3 Complete medium

Used to culture *P. falciparum* isolates or clones.

- I. Weigh 10.43g of powdered RPMI 1640 medium (RPMI 1640 with L-glutamine) into a 1L volumetric flask
- II. Add 5.95g HEPES (25mM final concentration)
- III. Add 42mL 5% sodium bicarbonate (as prepared in 6.1) (25mM final concentration)
- IV. Add:
  - a. 50mg hypoxanthine  
**Note:** Hypoxanthine should **not** be added if culturing *P.falciparum* for <sup>3</sup>H-hypoxanthine uptake inhibition assays ([INV07](#))
  - b. 430mg L-glutamine (0.3 g/L)
  - c. 1g glucose
  - d. 0.25mg orotic acid (0.25 mg/L; optional)
  - e. 10mg gentamycin or 32mg neomycin (10 mg/L; optional)
- V. Adjust volume to 1L with cell culture water
- VI. Stir until dissolved
- VII. Using a calibrated meter, check medium pH. If necessary, adjust to pH7.0 ± 0.3 with either hydrochloric acid (1N) or sodium hydroxide (1N). The medium should be orange, not red/purple or yellow, in colour due to phenol red dye
- VIII. Using a calibrated meter, check osmolarity = 292 ± 5% mOsm. If necessary, adjust with sterile water
- IX. In an aseptic environment, sterilize the medium using a 0.22µM filter
- X. Add sterile human serum to a final concentration of 10% or AlbuMAX™ I (lipid-enriched bovine serum albumin) to a final concentration of 0.5%.
- XI. Aseptically dispense medium into sterile bottles
- XII. Store at – 20 °C for a maximum 6 months
- XIII. After first use, store at 4 °C for one week
- XIV. Protect bottles from extended exposure to light

### 6.4 Quality control (QC)

Record:

- I. Preparation date
- II. Recorded powder weights
- III. Reagents' batch numbers
- IV. pH

QC records should be approved by a competent person and stored in a safe location.

## 7. References

Basco LK., *Field application of in vitro assays for the sensitivity of human malaria parasites to antimalarial drugs*. World Health Organisation. Available from: <http://www.who.int/malaria/publications/atoz/9789241595155/en/index.html> (Accessed: 29 November 2010)

Johnson JD, Denuall RA, GerenaL, Lopez-Sanchez M, Roncal NE, Waters NC. Assessment and continued validation of the malaria SYBR green I-based fluorescence assay for use in malaria drug screening. *Antimicrob Agents Chemoter* 2007 Jun; 51(6):1926-33. Abstract available from: <http://www.ncbi.nlm.nih.gov/pubmed/17371812> (Accessed: 29 November 2010).

Noedl H, Attlmayr B, Werbsdorfer WH, Kollaritsch H, Miller RS. A histidine-rich protein 2-based malaria drug sensitivity assay for field use. *Am J Trop Med Hyg* 2004 Dec; 71(6):711-4.

SOP Centre National de Référence du Paludisme, Paris, France.

SOP Institut de médecine tropicale du Service de santé des armées, Marseille, France.

SOP Institut Pasteur du Cambodge.

SOP Institut Pasteur de Guyane.

SOP Institut Pasteur de Madagascar